



Genesis™

Blood Collection Mixer

Key Features

- **Monitors collection volume to help reduce over/under bleeds (QNS units)**
- **Rocking action mixes blood with anticoagulant during collection**
- **Clamps tubing automatically at pre-set volume**
- **Operates on battery or A/C power**

The Genesis Blood Collection Mixer is an accurate, durable and reliable automated mixer that helps to reduce QNS rates and over bleeds in your operations by automatically clamping tubing when the preset volume amount is reached.

The Blood Collection Mixer operates off of a rechargeable battery or A/C power and is perfect for mobile and fixed sites.

Genesis™

Blood Collection Mixer

The Genesis Blood Collection Mixer monitors the amount of blood collected during a donation and automatically clamps tubing when the preset volume amount is reached. The mixer also features a rocking action that mixes blood with anticoagulant during the collection.



Key Features

- **Enhanced digital display provides key information: collected volume in mL, total weight in grams, elapsed time of donation**
- **Audible and visual alarms alert the end of collection, slow bleed, low battery**
- **Deep tray with tubing guides-holds all bags and in-line filters**
- **Magnetic tray for easy set-up and stability**
- **Operates on battery or A/C power**
- **Durable aluminum housing**
- **Tares Collection bag and filters**

Included items

- **1 year manufacturer warranty**
- **Heavy duty nylon carrying case**
- **Durable aluminum stand**
- **Rechargeable Ni-Cad battery and charger**
- **Power cord**
- **Operators manual**

Specifications

- **External dimensions: 11.5”L x 8.0” W x 5.5”H**
- **Collection Volume: 0 - 1199 mL**
- **Weight accuracy: +/- 1 gram**
- **Weight (mixer with battery): 7 lbs**

Optional accessories

- **Battery pack (part 425BN122C)**
- **6 unit multi-charger (part 425BC162)**

GENESIS Blood Collection Mixer is a trademark of Genesis BPS LLC
FENWAL is a trademark of Fenwal, Inc.

The Fenwal team of blood specialists are ready to help you make the most of your operations. Contact your local Fenwal Account Executive today or visit us at www.fenwalinc.com.
