



Genesis™

Blood Collection Monitor

Key Features

- **Monitors collection volume to help reduce over/under bleeds (QNS units)**
- **Clamps tubing automatically at pre-set volume target**
- **Operates on battery or A/C power**
- **Programmable default settings**

The Genesis Blood Collection Monitor electronically monitors and measures blood flow and volume during collection to help reduce QNS rates and over bleeds in your operations. The monitor automatically clamps the tubing when the preset volume amount is reached.

The Blood Collection Monitor operates off of a rechargeable battery or A/C power and is perfect for mobile and fixed sites.

Genesis™

Blood Collection Monitor

The Genesis Blood Collection Monitor electronically monitors and measures blood flow and volume during collection and automatically clamps the tubing when the preset volume amount is reached. The audible and visual alarms will alert the end of collection, slow bleed, low flow and provide a mix reminder. These features will help reduce costly over and under bleeds.



Key Features

- **Enhanced digital display provides key information: collected volume in mL, total weight in grams, elapsed time of donation**
- **Audible and visual alarms alert the end of collection, slow bleed, low flow, mix reminder**
- **Operates on battery or A/C power**
- **Performs automatic weight calibration**
- **Tares collection bag**

Included items

- **1 year manufacturer warranty**
- **Pole mount**
- **Rechargeable NiMH battery and charger**
- **Power cord**
- **500g standard calibration weight**
- **Operators manual**

Specifications

- **External dimensions: 6.4”L x 4.8” H x 2.4”W**
- **Collection Volume: 0 – 1199 mL range**
- **Weight accuracy: +/- 1 gram**
- **Weight (collection monitor with pole mount and a/c adaptor): 2.25 lbs**

Optional accessories

- **Heavy Duty aluminum case (part 425205013)**
- **Stand (pole and pedestal, part 425101052)**

GENESIS Blood Collection Monitor is a trademark of Genesis BPS LLC
FENWAL is a trademark of Fenwal, Inc.

The Fenwal team of blood specialists are ready to help you make the most of your operations. Contact your local Fenwal Account Executive today or visit us at www.fenwalinc.com.
