



Genesis™

Tube Sealers

Key Features

- **“Intelligent” RF Board - accepts and adjusts to all sizes of blood tubing**
- **Notched Seal™ for easy separation**
- **No warm-up required - instant on**
- **Durable aluminum housing**

The Genesis line of Tube Sealers are the industry standard for quality and versatility. The extensive line of dependable Genesis sealers provide options for bench-top units, portable units, space saving units and flexible multi-head solutions.

Find out more about the sealer line trusted by blood banks around the world.

Genesis™

Tube Sealers

The Genesis line of Tube Sealers are the industry standard for quality and versatility. The extensive line of dependable Genesis sealers provide options for bench-top units, portable units, space saving units and flexible multi-head solutions.



Name Part #	SE330 428SE330	SE340 428SE340	SE370 428SE370	SE530 428SE530	SE540 428SE540	SE640 428SE640
Description	Standard Benchtop Sealer	Multiple Station Benchtop Sealer	Battery Powered Benchtop Sealer	Space Saving Benchtop Sealer	Space Saving Benchtop Sealer	Portable Sealing Unit
Sealing Head	Fixed sealing head	Fixed & hand-held sealing head unit	Fixed & hand-held sealing head unit	Fixed sealing head	Fixed & hand-held sealing head unit	Hand-held sealing unit
Sealing time	0.5 – 1.2 seconds	0.5 – 1.2 seconds	0.5 – 1.2 seconds	0.5 – 1.2 seconds	0.5 – 1.2 seconds	0.5 – 1.2 seconds
“Intelligent” RF Board	Yes	Yes	Yes	Yes	Yes	Yes
Power	A/C	A/C	Battery Powered	A/C	A/C	Battery Powered
Splash guard	Retractable splash guard	Retractable splash guard	Retractable splash guard	Splash guard	Splash guard	N/A
Easy to clean head protector	Yes	Yes	Yes	Yes	Yes	Yes
Unit Weight	11.0 lbs	11.0 lbs	11.0 lbs	9.7 lbs	9.7 lbs	4.0 lbs
Special features	Durable aluminum housing	<ul style="list-style-type: none"> • 6' standard cable • Use fixed and hand –held sealing unit independently 	<ul style="list-style-type: none"> • 6' standard cable • Use fixed and hand –held sealing unit independently • 2000+ seals per full battery charge • Aluminum carrying case included 	<ul style="list-style-type: none"> • Retractable multi-head option • Custom sealer stand • Easy clean head protector • Slim design to maximize space • Flexible head spacing (75/80/85/90 mm) 	<ul style="list-style-type: none"> • Flexible multi-head option • 6' standard cable • Custom sealer stand • Slim design to maximize space • Flexible head spacing (75/80/85/90 mm) 	<ul style="list-style-type: none"> • 6' standard cable • Hand sealing unit • Battery operated and rechargeable • Lightweight canvas carrying bag included • 1000+ seals per full charge

GENESIS and NOTCHED SEAL are trademarks of Genesis BPS LLC
FENWAL is a trademark of Fenwal, Inc.

The Fenwal team of blood specialists are ready to help you make the most of your operations. Contact your local Fenwal Account Executive today or visit us at www.fenwalinc.com.