

Transfer Pack™ Containers




Transfer Pack containers in multiple configurations to address your needs




4R2001

4R2001 **Transfer Pack container with Coupler - 150 mL (bulk pack)**



4R2004


4R2004 **Quadruple Transfer Pack container with Coupler – 150 mL**



4R2023

4R2014 **Transfer Pack container with Coupler - 300 mL (bulk pack)**

4R2023 **Transfer Pack container with Coupler - 600 mL**



4R2027

4R2024 **Transfer Pack container with Male Luer – 600 mL**
Commonly used for therapeutic donations, empty bag without anticoagulant creates cost effective option for therapeutic collections, luer lock allows ease of connection to a needle

4R2027 **Transfer Pack container with 8 Couplers – 600 mL**
Eight couplers allow for pooling of products from multiple units

4R2255P **Quadruple Transfer Pack container Intended for Plasma Collection**
Quadruple Transfer Pack container intended for plasma collection with male luer adapter, one 1000 mL and three 300 mL containers

4R2335H **Product Sample Pack**
One plastic sample container for sterile connection with Amicus™ apheresis kits. Holds up to 5 mL.

Transfer-Pack, Amicus and Fenwal are trademarks of Fenwal, Inc.

The Fenwal team of blood specialists are ready to help you make the most of your operations. Contact your local Fenwal Account Executive today or call us at **1-800-3-Fenwal**.