

Transfusion Sets



Multiple options and configurations to address your transfusion needs



4C2160 Fenwal Blood Component Recipient Set with Standard Blood Filter and Luer Adapter

Standard 170-260 micron filter with female luer adapter for syringe, transfuses all components, approved for use up to 4 hours, 55" in length



4C2161 Fenwal Y-Type Blood Component Recipient Set with Standard Blood Filter and Luer Adapter

Standard 170-260 micron filter with Y-set, transfuses all components, approved for use up to 4 hours, 65" in length



4C2223 Blood Component Infusion Set

Blood component infusion set with female luer connector for syringe, transfuses all components, distal needle/catheter adapter contains 170-260 micron blood component filter, 13" in length

4C2240 2 Units Fenwal Plasma Transfer Set

Plasma transfer set with coupler and male luer lock, can connect to an IV or syringe with luer lock



4C2243 2 Units Fenwal Plasma Transfer Set

Plasma transfer set with two couplers.

4C2498H Y -Type Connector Set with Male Luer Lock

Commonly used for autologous transfusions with 4R0499



4R4402 Pressure Infusor, 600 mL

Provides constant pressure during transfusion

4R4403 Pressure Infusor, 1000 mL

Provides constant pressure during transfusion

FENWAL is a trademark of Fenwal, Inc.

The Fenwal team of blood specialists are ready to help you make the most of your operations. Contact your local Fenwal Account Executive today or call us at **1-800-3-Fenwal**.